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
EN Operation Manual

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Vacuson 40/60

Medical Suction Pumps

Vorwort

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Foreword

Congratulations on your purchase of a NOUVAG AG product. Thank you for the confidence shown in our products. Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years. You will find the conformity statement and list of authorized service representatives attached.

- **Please read instructions carefully before operating!**

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Table of contents

1	Product description	2
1.1	Intended use and operation	2
1.2	Contraindications	2
1.3	Technical data, Vacuson 40/60	2
1.4	Ambient conditions	2
1.5	Warranty coverage	2
2	Explanation of symbols	3
3	Safety information	3
3.1	EMC Manufacturer's Declaration of Conformity	3
3.2	Modification and misuse	3
3.3	Essential requirements	4
3.4	During use	4
4	Scope of delivery	4
5	Device overview	5
6	Startup	6
6.1	Device setup	6
6.2	Connection to the power supply	6
6.3	Preparation of secretion jars	6
6.4	Device preparation	7
7	Operation	9
7.1	Switching device on and off	9
7.2	Pump activation by On/Off-AIR-Pedal	9
7.3	Variegate with Vario-AIR-Pedal	9
7.4	Regulating suction process	9
7.5	Emptying secretion jar	10
7.6	Function control	10
8	Cleaning, disinfection and sterilization	11
8.1	Control unit and On/Off-AIR-Pedal	11
8.2	Secretion jar and jar lid	11
8.3	Bacteria filter	11
8.4	Silicone tubes	11
8.5	Cannulas and cannula handlebar	12
8.6	Quiver	12
9	Maintenance	13
9.1	Replacing the control unit fuse	13
9.2	Safety inspections	13
9.3	Bacteria filter	13
9.4	Secretion jar	14
9.5	Function control of float gauge valve	14
9.6	Disassembly of the Overflow Protection System	14
9.7	Seals and tubes	15
10	Malfunction and troubleshooting	16
11	Accessories and spare parts	17
12	Information on disposal	18

1 Product description

1.1 Intended use and operation

The Vacuson pumps are used for the suction of fluids and secretions. The pumps performance can be continuously regulated and monitored by the manometer.

The Vacuson suction pumps can be used in the following areas:

- Dental laboratory
- Medicine
- Intensive care
- Surgery
- Postoperative care

1.2 Contraindications

Relative or absolute contraindications may result from the general patient status or the risks of anesthesia related to the treatment.

1.3 Technical data, Vacuson 40/60

	Vacuson 40	Vacuson 60
Voltage:	115 V~at 60 Hz; 230 V~at 50 Hz	115 V~at 60 Hz; 230 V~at 50 Hz
Power consumption:	230 VA	320 VA
Fuses:		
for 115 V model:	2 x T2 A, L, 250 VAC	2 x T2 A, L, 250 VAC
for 230 V model:	2 x T1 A, L, 250 VAC	2 x T1 A, L, 250 VAC
Protection class:	Class I	Class I
Applied part:	Type BF	Type BF
Adjustable vacuum:	– 0.9 bar at 686 mmHG	– 0.9 bar at 686 mmHG
Dimensions, W x H x D:	360 x 300 x 280 mm	360 x 300 x 280 mm
Weight:	10 kg	14 kg
Accuracy limit, Manometer:	5 %	5 %
Suction Pump Capacity:	40 l/min (for watery solutions, without attached instrument)	60 l/min (for watery solutions, without attached instrument)
Pump motor power:	185 Watt	250 Watt

1.4 Ambient conditions

	Transport and storage:	Operation:
Relative humidity:	Max. 90 %	Max. 80 %
Temperature:	0°C – 60°C	10°C – 30°C
Atmospheric pressure:	700 hPa – 1060 hPa	800 hPa – 1060 hPa

1.5 Warranty coverage

Purchasing a Vacuson 40 or Vacuson 60 suction pump entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further **6 month**.

Consumable parts are not covered by the warranty. Improper use or repair, or failure to observe these instructions, relieve us from any obligations arising from warranty provisions or other claims.

2 Explanation of symbols

	Important information		Autoclavable at max. 135°C for at least 5 min.
	Warning		Suitable for thermal disinfection
	Date of manufacturing		Protective ground
	Manufacturer		Observe the instructions for use
	Type BF applied part is the filling tube with connected instruments		Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Valid local disposal regulations apply.
IPX8	Protection continual submerging.	SN	Symbol indicating the serial number with the date of manufacture (year/month).
	Pedal	REF	Symbol indicating the order number.
	Date of expiry	LOT	Symbol indicating the lot number.
CE 1275	CE symbol with notified body		Not for reuse

EN

3 Safety information

Your safety, the safety of your team, and of course that of your patients is very important to us. It is therefore essential to bear the following information in mind:

3.1 EMC Manufacturer's Declaration of Conformity

Please observe the information on electromagnetic compatibility provided with these operating instructions.

3.2 Modification and misuse

- Modification or manipulation of the Vacuson 40/60 suction pumps and its accessories is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be canceled.
- Use of the Vacuson 40/60 suction pumps outside the indications described in Section 1.1 is prohibited. The user or operator is solely responsible for any such use.

3.3 Essential requirements

The Vacuson 40/60 suction pumps may only be operated under constant supervision of qualified and trained personnel!

Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims!

The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.

Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.

Repairs may only be performed by authorized NOUVAG service technicians.

3.4 During use

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (refer to chapter 8 "Cleaning, disinfection and sterilization").

At the choosing of the instrument the operator has to make sure it's biocompatible, according to EN ISO 10993.

While in operation the control unit of the Vacuson 40/60 suction pumps must be at least 1 meter above ground.


Do not use the device in the vicinity of flammable mixtures.

4 Scope of delivery

Ref.	Description	Quantity
	Vacuson 40 (Ref. 4227-115 V/4227-230 V)	
4275	Control unit Vacuson 40 incl. On/Off-pneumatic-Pedal -----	1
4075	Connection tube 8 x 3 x 400 mm between bacteria filter and jar lid, silicone, sterilizable -----	1
4076	Filling tube 8 x 3 x 1700 mm, silicone, sterilizable including bacteria filter -----	1
	Vacuson 60 (Ref. 4237-115 V/4237-230 V)	
4280	Control unit Vacuson 60 incl. On/Off-pneumatic-Pedal-----	1
4075	Connection tube 8 x 3 x 400 mm between bacteria filter and jar lid, silicone, sterilizable -----	1
4076	Filling tube 8 x 3 x 1700 mm, silicone, sterilizable including bacteria filter-----	1
	Optional:	
4122	Stand for Vacuson -----	1
4242	Vario-AIR-Pedal -----	1
4052	Secretion jar, 2 liter, polypropylene -----	1
4245	Secretion jar, 5 liter, polypropylene -----	1
4058	Jar lid for 2 and 5 liter secretion jars-----	1
4035	2 liter disposable secretion pouches including lid for MONOKIT system -----	50
4036	2 liter polycarbonate receptacle for MONOKIT system-----	1
4037	Bottle holder for MONOKIT system for mounting secretion jars on Vacuson pump-----	1
31997	Operating instructions on CD-ROM -----	1


5 Device overview

Front view

1. Cannula handlebar (optional)
2. Cannula (optional)
3. Quiver with mounting bracket (optional)
4. Filling tube (170 cm), silicone
5. Stand with antistatic rolls and brakes
6. On/Off pneumatic pedal
7. Connection tube of pneumatic pedal, one piece with pedal
8. Secretion jar (Example 2 Liter)
9. Jar lid of secretion jar
10. Secretion jar mount
11. Connection tube (40 cm), Silicone
12. Intake nozzle with bacteria filter
13. Carrying handle
14. Manometer
15. Ready indication, LED
16. Vacuum regulator (VACUUM)
17. Mount for Quiver or secretion jar
18. Type plate with type designation, reference number, serial number, information on power supply and device fuse.
19. Port for pneumatic pedal, On/Off
20. Air exhaust port (EXHAUST)
21. Plug for VARIO-AIR pedal
22. Power plug socket
23. Main switch
24. Fuse compartment

Rear view

6 Startup

6.1 Device setup

- Place the Vacuson 40/60 suction pump and all required accessories and instruments on an even, non-slip surface and make sure you have good access to all controls.
- Do not allow the operating range of the device (including cable) and the connected instrument to be compromised by limiting factors.
- The manometer must be fully visible at all times.
- The On/Off-AIR-Pedal must be placed within stepping distance between the patient and the surgeon.
- It must be explicitly ensured that no objects can fall on the pedal.
- The power plug at the rear of the device must be accessible at all times.
- The ventilation slots at the housings bottom and sideways of the Vacuson 40/60 must be kept clear in order to prevent temperature from becoming excessive.
- While in operation the Vacuson 40/60 suction pumps must be at least 1 meter above ground. Nouvag recommends for convenient handling of the device the use of the optional movable stand.

6.2 Connection to the power supply

Before switching on, make sure that the power supply unit of the device matches the country's specific service voltage!

The power supply unit of the Vacuson 40/60 pumps is not swichable to the country specific service voltage. The device has to be ordered according to the country specific service voltage.

In order to prevent the risk of an electric shock, the device may only be connected to a power network with a PE protective ground conductor.

The power plug socket is located at the rear of the device.

6.3 Preparation of secretion jars

1. Hold open jar available (2 or 5 liter).
2. Press Jar Lid with Turn and Tilting Lever in open-position firmly onto the jar (the latch of the locking system is in open-position).

3. Rotate Turn and Tilting Lever by 180° (Turn and Tilting Lever now facing away from the grasp). Make sure the gripper catches the rim of the jar.

4. Flap down Turn and Tilting Lever into the designated groove.


Turn and Tilting Lever locked and secured.

EN

6.4 Device preparation

1. Prior to use all sterilizable parts (Tubes, Cannulas, Cannula handlebar, char lid and bottles) must be sterilized.

2. Attach secretion jars with mounted and locked lid to the device.


3. Lay the pneumatic pedal on the floor and plug the connection tube into the port for the pneumatic pedal.

4. The Vario-AIR-Pedal is connected to the air inlet port at the control unit by the connection tube. If not in use the port is covered with a lid. Make sure to cover the port again when Vario-AIR-Pedal is not connected.

5. Attach quiver (optional) at the mounts of the rollable stand (optional).

6. Attach the short connection tube (40 cm) on one end with the bacteria filter and the other end with the narrow connector (VACUUM) of the jar lid.

7. Attach the filling tube on one end with the angled connector (optional Ref. 28535) and the other end of the tube with the instrument.

8. Mount connection tube with bacteria filter onto the intake nozzle at the top of the Vacuson pump. Graft the other end of the connection tube with the narrow connector onto the smaller nozzle of the secretion jar (VACUUM).

9. Graft the angled connector (optional Ref. 28535) of the filling tube (170 cm) onto the wider nozzle (PATIENT) of the secretion jar. Hang the other end of the filling tube with the instrument into the quiver.

10. Screw suction cannula (optional) onto the cannula handlebar (Vacuson 60, Liposuction) and hang it back into the quiver.

11. Connect the device socket with the wall socket, using the device's power cord.

Before switching on, make sure that the power supply of the device matches the country's specific service voltage!

7 Operation

7.1 Switching device on and off

Use the power switch "I/O", at the rear of the device, to switch the device on and off. The standby is signaled by the LED status light at front of the device. The device can be switched off at any time, irrespective of any procedure for device switch-off.

EN

7.2 Pump activation by pneumatic pedal

The included standard pneumatic pedal when pushed activates a switch in the Vacuson pump, to switch the pump on and off. This is maintained by an air pillow. After the device was switched on, the pneumatic pedal has to be pressed to activate the pump and generate vacuum. Pressed again the pump is deactivated. The pump can only be activated by the pedal. The last used state before switch-off stays active.

7.3 Variegate with Vario-AIR-Pedal

Handsfree regulation of the vacuum can be maintained by using the Vario-AIR-Pedal (optional, Art-No. 4242). When pressed, the Vario-AIR-Pedal opens a valve and the vacuum is reduced by forced ventilation of the pressure system. The more the pedal is pressed the less vacuum can build up, hence the suction performance decreases. If the suction performance is controlled by the Vario-AIR-Pedal, the Vacuum Controller on the front side of the device ideally remains in maximal position. The Vario-AIR-Pedal is normally used together with the pneumatic pedal.

7.4 Regulating suction process

The suction process is regulated by the Vacuum controller at the front side of the device.

<i>Rotate clockwise:</i>	<i>Vacuum increases, suction performance is enhanced.</i>
<i>Rotating anticlockwise:</i>	<i>Vacuum decreases, suction performance is reduced.</i>

Establishing airtightness of the suction system:

1. Turn Vacuum Controller anticlockwise all the way to the stop (equivalent to vacuum minimum).
2. Switch on mains switch (I/O) of the pump and press On/Off pedal shortly. Pump is running and building up vacuum.
3. Crimp suction tube to generate maximal air-tightness of the suction system.
4. Turn Vacuum Controller clockwise (equivalent to vacuum maximum).
5. Wait for maximum build up of vacuum (equivalent to ≥ -0.9 bar).
6. Now by turning the Vacuum Controller the suction performance can be regulated steplessly from 0 to -0.9 bar.

7.5 Emptying secretion jar

The jar lid of the secretion jar is equipped with an overflow protection system to prevent the vacuum system from being flooded by secretion fluids at high filling levels of the secretion jar. Therefore a float gauge is responsible. At high filling levels of the secretion jar and the resulting locking of the overflow protection system the secretion jar has to be emptied or replaced by another secretion jar.

1. Switch off suction pump.
2. Disconnect tubes from the secretion jar lid.
3. Unhinge full secretion jar from the secretion jar mount and dispose of secretion fluids according to national disposal regulations.
4. Procure used secretion jars to the reprocessing cycle.
5. Connect suspended tubes with new, ready to use secretion jar.

7.6 Function control

To obtain trouble-free operation of the suction pump all the components and functions of the pump system have to be tested prior every assignment.

General functions:

1. Use the power switch "I/O", at the rear of the device, to switch the device on, LED is illuminated.
2. Device fan is running.

Suction pump with pneumatic pedal:

1. Use the power switch "I/O", at the rear of the device, to switch the device on, LED is illuminated.
2. Device fan is running.
3. Hold hand in front of the «Exhaust» (Rear of the device). Airflow is perceptible.
4. Press pneumatic pedal shortly to activate the suction pump.
5. Turn Vacuum Controller clockwise (equivalent to vacuum maximum).
6. Control vacuum intensity at the opening of the cannula. Strong suction performance.
7. Turn Vacuum Controller anticlockwise all the way to the stop (equivalent to vacuum minimum).
8. Control vacuum intensity at the opening of the cannula. Weak suction performance.
9. Use the power switch "I/O", at the rear of the device, to switch off device, LED is not illuminated.

Suction pump with pneumatic pedal and Vario-AIR-Pedal:

1. Use the power switch "I/O", at the rear of the device, to switch the device on, LED is illuminated.
2. Device fan is running.
3. Hold hand in front of the «Exhaust» (Rear of the device). Airflow is perceptible.
4. Press pneumatic pedal shortly to activate the suction pump.
5. Turn Vacuum Controller clockwise (equivalent to vacuum maximum).
6. Control vacuum intensity at the cannulas opening. Strong suction performance.
7. Press Vario-AIR-Pedal. The more it is pressed, the weaker the suction performance at the cannula.
8. Turn Vacuum Controller anticlockwise all the way to the stop (equivalent to vacuum minimum).
9. Control vacuum intensity at the cannulas opening. Weak suction performance.
10. Press Vario-AIR-Pedal. The more it is pressed, the weaker the suction performance at the cannula. The suction performance now is not perceptible anymore.
11. Use the power switch "I/O", at the rear of the device, to switch off device, LED is not illuminated.

Malfunctions and troubleshooting:

To solve problems refer to chapter 10 "Malfunctions and troubleshooting".

8 Cleaning, disinfection and sterilization

The following points in particular are important with regard to caring for the material:

- Perform cleaning, disinfection and sterilization after every treatment!

- Always autoclave the material in sterilization packaging.

- Make sure that sterilization packaging is no more than 80 % full.

- Always autoclave the material at max. 135°C for at least 5 minutes.

- If sterilized material is not used immediately, the material packaging must be labeled with the sterilization date.

- Novvag AG recommends including a sterility indicator.

EN

8.1 Control unit and pneumatic pedal

Control unit and pneumatic pedal do not come into contact with the patient.

Wipe the outside using micro-biologically tested surface disinfectant or a 70 % isopropyl solution. The front plate of the control unit is sealed for this purpose and can be wiped clean.

8.2 Secretion jar and jar lid

The reprocessing instructions for the secretion jar and jar lid is provided in the operating instructions delivered with the secretion jar.

8.3 Bacteria filter

The bacteria filter located on top of the Vacuson pump is a one way product and cannot be cleaned or sterilized.

A periodical replacement of the bacteria filter is recommended after 8 hours of use, but definitely after it came into contact with foam or infectious material.

After contact with watery solutions the bacteria filter locks down, because of its hydrophobic characteristic, to protect the pump from clogging. Hence the further operation of the pump is not possible. The bacteria filter has to be replaced.

8.4 Silicone tubes

Ref. 4075, connection tube 8 x 3 x 400 mm from bacteria filter to secretion jar, silicone, sterilizable

Ref. 4076, filling tube 8 x 3 x 1700 mm, from secretion jar to suction cannula, silicone, sterilizable

Reprocessing restrictions	Frequent reprocessing of the silicon tubes has only a limited impact. The end of the product service life is normally determined by wear and damage through use.
---------------------------	--

INSTRUCTIONS	
At location of use	No special requirements.
Storage and transport	No special requirements. Long holding times before reprocessing have to be avoided due to surface drying.
Preparation for cleaning	No special requirements.

Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures the connection of tubes to the washer-disinfector for rinsing. Use only neutral cleaning agents for this purpose. 1) Place silicone tubes in the load carrier. 2) Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with fully deionized water. 3) Perform a 10-minutes rinse cycle at 93°C to facilitate thermal disinfection. 4) When removing, check silicone tubes, to verify whether soiling is still visible. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: Neutral cleaning agent, soft brush, running, demineralized water (< 38°C) Procedure: 1) Rinse off and brush away surface soiling from the silicone tubes. 2) Rinse silicone tubes thoroughly under running water.
Manual disinfection	For manual disinfection, submerge silicone tubes in chlorinefree disinfection solution.
Drying	Allow silicone tubes to dry sufficiently in a drying cabinet.
Inspection and maintenance	Perform a visual inspection to check for damage, corrosion and wear.
Packaging	Individual: Pack silicone tubes in individual packaging for sterile items. Sets: Sort silicone tubes on trays intended for this purpose or place them on allpurpose sterilization trays.
Sterilization	 Autoclave in vacuum autoclave at max. 135°C for at least 5 minutes. When sterilizing several items during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in case of autoclaves without a post-vacuum function. Allow the silicone tubes to dry in the bag for at least one hour at room temperature with the paper side facing upwards. * Temperature exposure times are based on country-specific guidelines and standards.
Storage	No special requirements. If sterilized silicone tubes are not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.

The effectiveness of the sterilization instructions provided above for reprocessing this medical product has been validated by Nouvag AG. The user is responsible for ensuring that the sterilization procedure performed achieves the required results. This requires validation and routine monitoring of the procedure. The staff member who completes the procedure bears sole responsibility for any deviation on his part from the instructions provided. Deviations necessitate revalidation of the effectiveness of the procedure as well as of the technical resilience of the reprocessed items with regard to the modified sterilization process.

- The tube set Ref. 6024 (optional) is delivered in sterile condition. It is determined for single use and may not be resterilized!
- Contaminated tube sets have to be disposed of expertly!

8.5 Cannulas and cannula handlebar

The optional cannulas and the cannula handlebar are in contact with the patient and therefore have to be reprocessed adequately.

The reprocessing instructions are in the operation instructions, delivered together with the cannula and handlebar.

8.6 Quiver

Clean quiver from debris and soiling. Use a clean, damp cloth and/or an appropriate brush with disinfection agent.

1. Attention, it's important to use a disinfection agent compatible with polycarbonate.
2. Pack quiver in individual packaging for sterile items (siehe DIN 58953).
3. Autoclave wrapped quiver **at max. 135°C for at least 5 minutes***.
4. A drying cycle must be added in case of autoclaves without a post-vacuum function. Allow quiver to dry in the bag for at least one hour at room temperature with the paper side facing upwards.

If sterilized quiver is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.


** Temperature exposure times are based on country-specific guidelines and standards.*

9 Maintenance

9.1 Replacing the control unit fuse

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch:

- Unplug the power plug.
- Open the fuse slot using a screw driver.
- Replace defective fuses T 2 A, L, 250 VAC (115 V model)/ T 1 A, L, 250 VAC (230 V model).
- Slide the fuse holder back in and close the fuse slot.
- Plug in the power plug again.

1. Fuse slot locking mechanism
2. Fuse slot
3. Fuse 1
4. Fuse 2

9.2 Safety inspections

The performance of safety inspections on medical devices is required by law in several countries. The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of **2 years** applies to the Vacuson 40/60.

NOUVAG AG offers a safety inspection service for its customers. Addresses can be found in the appendix of this operation manual under “Service centers”. For further information please contact our technical service department.

Further international service centers are listed on the Nouvag website:

www.nouvag.com > Service > Service centers

9.3 Bacteria filter

A periodical replacement of the bacteria filter is recommended after 8 hours of use, but definitely after it came into contact with foam or infectious material. For reordering refer to chapter 11 to retrieve the article number.

9.4 Secretion jar

The influxing mixture of air and secretion fluids into the secretion jar causes the build up of foam. It's recommended to use an antifoam agent to suppress the build up of foam. Prior use of the secretion jar fill an anti foam agent into the clean, dry jar. Don't use disinfection solution, because most of them benefit the build up of foam.

Make sure the secretion jars are in good condition. Check the jars routinely for cracks and rifts and be sure the jars flange is immaculate. It's important to guarantee full air tightness of the system which is responsible for troublefree operation of the pump.

9.5 Function control of float gauge valve

The proper functioning of the overflow protection system, built in the jar lid, has to be checked periodically.

1. Connect jar lid (VACUUM) with bacteria filter, using the connection tube (8 x 3 x 400 mm).
2. Turn Vacuum Controller clockwise (equivalent to vacuum maximum) all the way to the stop.
3. Press pneumatic pedal to generate vacuum.
4. Press the float gauge of the overflow protection system towards the lid.
5. The manometer shows increasing values up to the maximum. (> -0.9 bar).

If the manometer doesn't show maximal vacuum (> -0.9 bar), the overflow protection system has to be disassembled, cleaned and the seals have to be replaced.


9.6 Disassembly of the Overflow Protection System

1. Unscrew threaded Overflow Protection connector.
2. Remove seal from inside connector (O-Ring acc. picture).
3. Clean Overflow Protection System and float gauge.
4. Install new seal (O-Ring acc. picture).
5. Reassemble of Overflow Protection connector.

Function control after reassembly:

6. Hold lid perpendicularly.
7. Press float gauge repeatedly towards the lid.
8. Float gauge must fall back in place by itself.

If the float gauge doesn't fall back in place by itself the cleaning procedure has to be repeated and the sitting of the O-Ring seal has to be checked and corrected.


9.7 Seals and tubes

To insure proper function of the Suction Pump, all tubes and seals must be periodically checked, and replaced after at most 250 sterilization cycles or five years of operation. Only when in perfect order can sufficient vacuum be built up.

Defective seals and tubes must be replaced.

EN

	1. Seal (O-Ring) at intake nozzle of control unit for connection with bacteria filter ----- 2 units ----- Ref. 4063
	2. Connecting tube between bacteria filter and secretion jar lid (VACUUM) ----- 1 unit ----- Ref. 4075
	3. Seal (O-Ring) at connection nozzle (unscrewable, VACUUM) at secretion jar lid ----- 4 units ----- Ref. 4064
	4. Seal (O-Ring) at overflow protection system at jar lid ----- 1 unit ----- Ref. 28958
	5. Main seal between jar lid and jar ----- 1 unit ----- Ref. 28957
	6. Seal (O-Ring) at connection nozzle (unscrewable, PATIENT) of secretion jar lid ----- 4 units ----- Ref. 4063
	7. Filling tube between connection nozzle (PATIENT) and cannula ----- 1 unit ----- Ref. 4076

10 Malfunction and troubleshooting

Malfunction	Cause	Solution	Reference in manual
Device is not operational	Pump is not switched on	Switch main switch "I/O" to "I"	7.1 Switching device on and off
	No connection to electricity-supply	Connect power cord to electricity-supply	6.4 Device preparation
	Wrong voltage	Check power supply of your Pump	5.0 Device overview
	Defective fuses	Replace fuse	9.1 Replacing the control unit fuse
Pedal is not functioning	On/Off-AIR-Pedal is not connected	Connect On/Off-AIR-Pedal with device at rear	6.4 Device preparation
	Control unit is not switched on	Switch main switch "I/O" to "I"	7.1 Switching device on and off
	Incorrect operation	Read instruction manual carefully	
Suction pump is not functioning	On/Off-AIR-Pedal is not connected	Connect On/Off-AIR-Pedal with device at rear	6.4 Device preparation
	Vacuum-system is not air tight	Check all seals and tubes. Make sure char lid is properly closed	6.3 Preparation of secretion jar 6.4 Device preparation 9.7 Seals and tubes
	Air inlet port at rear of the device is open	Close air inlet port with its cap.	6.4 Device preparation
	Tubes are connected wrong	Connect tubes correctly	6.4 Device preparation
	Jar is full and Overflow protection has locked down	Replace full jar by a fresh, empty jar	7.5 Emptying secretion jar
	Incorrect operation	Read instruction manual carefully	
Suction pump is not working properly	Vacuum controller is not opened wide enough	Turn Vacuum controller clockwise	7.4 Regulating suction process
	Vacuum-system is not air tight	Check all seals and tubes. Make sure char lid is properly closed	6.3 Preparation of secretion jar 6.4 Device preparation 9.7 Seals and tubes

If a fault cannot be rectified, please contact your supplier or an authorized service center. The addresses are provided on the last page of this operating instructions.

11 Accessories and spare parts

Seals

Ref.

(Refer to 9.7 Seals and tubes)

Connection elements

Ref.

	Standard straight wide connector (PATIENT) at secretion jar lid, attached to filling tube -----	4056
	Standard straight narrow connector (VACUUM) at secretion jar lid, attached to connecting tube -----	4047
	Angled connector (VACUUM) for more convenient tube connection with filling tube -----	28535

Tube elements

Ref.

	Connecting tube, Silicone, bacteria filter to jar lid, (Outer-Ø x inner-Ø x length) 8 x 3 x 400 mm -----	4075
	Filling tube, Silicone, Autoclavable (Outer-Ø x inner-Ø x length) 8 x 3 x 1700 mm -----	4076
	Single use filling tube Polypropylene, sterile, length 4 m -----	6026

Accessories

Ref.

	Secretion jar, Polypropylene, unbreakable, autoclavable, 2 liter volume -----	4052
	Secretion jar, Polypropylene, unbreakable, autoclavable, 5 liter volume -----	4245
	MONOKIT, disposable secretion jar system with receptacle for 2 liter pouches, with suspension device -----	4038
	MONOKIT, 2 liter disposable pouches including lid, 50 units -----	4035
	Bacteria filter for all NOUVAG Vacuson suction pumps, disposable, hydrophobe -----	4246
	Secretion jar lid with overflow protection system for 2 and 5 liter secretion jars, complete -----	4058
	Quiver, sterilizable, 30 cm length, with suspension device -----	4043
	Quiver, sterilizable, 40 cm length, with suspension device -----	4044
	Yankauer suction cannula, length 28 cm, Ø 2.0 mm -----	4446
	Andrews cannula, length 24 cm, Ø 2.0 mm -----	4449
	Suction cups for obstetrics with Ø 60 mm -----	4053
	Suction cups for obstetrics with Ø 70 mm -----	4054
	Tubing holder to mount on SANICAR or STATIV -----	4118

Suction Cannulas for Liposuction

Ref.

	Cannula handlebar with opening for false air ventilation -----	4391
	Cannula handlebar without opening -----	4390
	Curved cannula for Femoral Liposuction, Ø 3 mm, length 200 mm, 22 openings 1.5 mm -----	4362
	Curved cannula for Femoral Liposuction, Ø 3 mm, length 300 mm, 30 openings 1.5 mm -----	4365
	Curved cannula for Femoral Liposuction, Ø 4 mm, length 200 mm, 22 openings 1.5 mm -----	4368
	Curved cannula for Femoral Liposuction, Ø 4 mm, length 300 mm, 30 openings 1.5 mm -----	4372
	Angled cannula, 30°, for Femoral Liposuction, Ø 3 mm, length 200 mm, 22 openings 1.5 mm -----	4381
	Straight cannula, Ø 1.5 mm, length 150 mm, 1 oval opening -----	4361
	Straight cannula, Ø 2 mm, length 150 mm, 1 oval opening -----	4364
	Straight cannula, Ø 2 mm, length 150 mm, 18 opening 1 mm -----	4373
	Straight cannula, Ø 3 mm, length 150 mm, 18 openings 1.5 mm -----	4374
	Straight cannula, Ø 3 mm, length 200 mm, 22 openings 1.5 mm -----	4378
	Straight cannula, Ø 3 mm, length 300 mm, 30 openings 1.5 mm -----	4387
	Straight cannula, Ø 4 mm, length 200 mm, 22 openings 2.0 mm -----	4379
	Straight cannula, Ø 4 mm, length 300 mm, 30 openings 1.5 mm -----	4388

Infiltration cannula

Ref.

	Straight cannula, Ø 3 mm, length 250 mm -----	4350
---	---	------

To order additional parts, please contact our customer service department.

12 Information on disposal

When disposing of the device, device parts and accessories, the regulations prescribed by law must be observed.

Do not dispose of devices with household waste!

To ensure environmental protection, old devices can be returned to the dealer or manufacturer.

Motors that have reached the end of their service life may not be disposed of with household waste. Motors must be sterilized before disposal. Please observe currently valid national disposal regulations for infectious waste.

Contaminated single-use tubing sets are subject to specific disposal requirements. Please observe currently valid national disposal regulations for infectious waste.

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Electromagnetic compatibility (EMC)

Remark:

The **Product** subsequently referred to herein always denotes the Vacuson 40 / 60.

Changes or modifications to this product not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product shall not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the product shall be tested to verify normal operation in the configuration in which it is being used.

Compliant Cables and Accessories

WARNING

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product.

The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Description	Length max.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Product requires continued operation during power mains interruptions, it is recommended that the Product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	Not applicable	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity for not life support equipment			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V rms 150 kHz to 80 MHz outside ISM bands	3 V rms 150 kHz to 80 MHz outside ISM bands	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			where P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.		
b	over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment			
The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

KONFORMITÄTSERKLÄRUNG / DECLARATION OF CONFORMITY / DECLARATION DE CONFORMITE / DICHIARAZIONE DI CONFORMITA / DECLARACIÓN DE CONFORMIDAD

Wir, die Firma
We, of the company
Nous, la firme
Noi, della ditta
Nosotros, la empresa

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
For all Nouvag servicecenters worldwide please check:

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LGA InterCert GmbH
Nouvag AG • St.Gallerstr. 23-25 • CH-9403 Goldach
Nouvag GmbH • Schulthaißstrasse 15 • D-78462 Konstanz
Nouvag USA Inc. • 18058 Albyn Court • Lake Hughes, CA 93532 • USA